



CR-TECH DICOM Conformance Statement

Last update of this document: 02/01/07

Introduction

This section is an abbreviated DICOM conformance statement for CR-TECH products.

Supported transfer syntaxes (Reading)

The Transfer Syntax UID is in the file DICOM Tag field (0002,0010).

	Description
Uncompressed Transfer Syntax	
1.2.840.10008.1.2	Implicit VR - Little Endian
1.2.840.10008.1.2	Modality Worklist
1.2.840.10008.1.2.1	Explicit VR - Little Endian
1.2.840.10008.1.2.1.99	Deflated Explicit VR - Little Endian
1.2.840.10008.1.2.2	Explicit VR - Big Endian
RLE Transfer Syntax	
1.2.840.10008.1.2.5	Description Run Length Encoding, Lossless
JPEG Transfer Syntax	
1.2.840.10008.1.2.4.50	Description JPEG Baseline (Process 1)
1.2.840.10008.1.2.4.51	JPEG Extended (Process 2 & 4)
1.2.840.10008.1.2.4.57	JPEG Lossless, Non-Hierarchical (Process 14)
1.2.840.10008.1.2.4.70	JPEG Lossless, Hierarchical, First-Order Prediction (Process 14, [Selection Value 1])

Supported transfer syntaxes (Writing)

The Transfer Syntax UID is in the file DICOM Tag field (0002,0010).

Uncompressed Transfer Syntax	Description
1.2.840.10008.1.2	Implicit VR - Little Endian
1.2.840.10008.1.2.1	Explicit VR - Little Endian
JPEG Transfer Syntax	
1.2.840.10008.1.2.4.70	Description JPEG Lossless, Hierarchical, First-Order Prediction (Process 14, [Selection Value 1])



Supported "Photometric Interpretation" pixel format (Reading)

The Photometric Interpretation UID is in the file DICOM Tag field (0028,0004).

Photometric Interpretation" pixel format	Description
MONOCHROME1	grey level image description
MONOCHROME2	grey level image description
PALETTE COLOR	pseudo color image description
RGB	true color image description
YBR_FULL	true color image description
YBR_FULL_422	true color image description

Supported "Photometric Interpretation" pixel format (Writing)

The Photometric Interpretation UID is in the file DICOM Tag field (0028,0004).

Photometric Interpretation" pixel format	Description
MONOCHROME2	grey level image description (high values=bright, low values=dark)

Supported 'Bits Allocated' values (Reading)

The Bits Allocated value is in the file DICOM Tag field (0020,0100).

Classical values	Description
8, 10,12,14,15, 16	12 means that 1 pixels is stored in 1 'short int' or 2 bytes



Supported SOP Classes - User Of Service (SCU) and Provider Of Service (SCP)

Networking SOP Classes	Description
Verification	SCU & SCP
Image Transfer	SCU & SCP, all modalities, multicast
Query & Retrieve	Study Root Q/R - FIND, Study Root Q/R - MOVE, SCU & SCP, multicast
Modality Worklist	Information Model SCU

Support of Extended Character Sets

Configurable Parameters:

The following items are configurable for each entity in the configuration file:

- Port Number
- AE Title for SCU operations
- AE Title for Storage SCP
- AE Title for Query/Retrieve SCP / SCU

- AE Title for Modality Worklist



Processing Application Support

SOP Classes

Description



Digital X-Ray Image Storage - for Presentation
Digital X-Ray Image Storage - for Presentation
Digital X-Ray Image Storage - for Processing Storage
Digital Mammography X-ray Image Storage - For Presentation
Digital Mammography X-ray Image Storage - For Processing
Digital Intra-Oral X-ray Image Storage - For Presentation
Digital Intra-Oral X-Ray Image Storage - For Processing
CT Image Storage
US Multi-frame Image Storage
MR Image Storage
Enhanced MR Image Storage
MR Spectroscopy Storage
Ultrasound Image Storage
Secondary Capture Image Storage
Multi-Frame Single Bit Secondary Capture Image Storage
Multi-Frame GrayScale Byte Secondary Capture Image Storage
Multi-Frame GrayScale Word Secondary Capture Image Storage
Multi-Frame True Color Secondary Capture Image Storage
Standalone Overlay Storage
Standalone Curve Storage
Standalone Modality Worklist Storage
Standalone VOI LUT Storage
GrayScale Softcopy Presentation State Storage SOP Class
X-Ray Angiographic Image Storage
X-Ray Radiofluoroscopic Image Storage
Nuclear Medicine Image Storage
Raw Data Storage
Positron Emission Tomography Image Storage
Standalone PET Curve Storage

Use and Create

Print SOP Classes - User Of Service (SCU)

SOP Classes

SOP Class UID



Print Job SOP Class

1.2.840.10008.5.1.1.14

Printer SOP

1.2.840.10008.5.1.1.16